



210

ACTIS BROACH

ACCESSORY

Length: 97mm
Author: Robert Judet

EST.
1947



Accessory Actis Broach #0

Common Uses:

Hip replacement surgeries

Primary Properties:

- Extraction Teeth
- Compaction Teeth
- Distal Cutting Teeth
- Hydroxyapatite layered stem

Function:

The Accessory Actis Broach enables many surgical approaches for hip replacement which allows for the preservation of patient anatomy while ensuring adequate fit in the canal.

Variations:

- Titanium sintered beads
- Collarless/Collared
- Lengths: 97mm - 119mm (#0 - #11)

Did you know?



The implant and instrumentation are designed to facilitate proper fit and placement of the femoral prosthesis using tissue sparing techniques.

